

Savina Technical Service Bulletin # 2

Re:	New software version 1.11
Date:	July 31, 2001
Reference Doc:	IDM Savina # 2, dated May 08, 2001
Reason:	Dräger Medizintechnik has released software version 1.11.
Devices affected:	The software update with software version 1.11 is mandatory for all Savina units with installed software versions 1.00, 1.01 and 1.10.
When:	ASAP, prior to shipment to Customer
Reporting:	To Dräger Medical, Inc. by Fax 1-215-721-5789. A Feedback Form is included in the conversion description.
Additional Info :	<p>For Savina with installed SW version 1.00 only:</p> <ul style="list-style-type: none">• First install SW 1.01 (part of the SW 1.11 kit)• Power on unit with SW 1.01 and generate a device ID number by changing all fields for date and clock in the configuration menu 3.• If the ID number in the configuration menu 4 is available set date and time back to the correct values.• Install SW 1.11.
Cost:	N/A
Ordering Info:	SW 1.11 Kit Savina P/N 84 14 535
Distribution:	DrägerService Personnel and authorized Service Organizations for CCS products.

If you have any questions, please contact Technical Support by phone at 1-800-543-5047 or by fax at 1-215-721-5789.

Dräger Medical, Inc.
Technical Product Manager